

# Prodigy™ Phenomenex

充てん剤	粒子径 (μm)	細孔径 (Å)	細孔容量 (mL/g)	表面積 (m <sup>2</sup> /g)	炭素含有率 (%)	表面被覆率 (μmole/m <sup>2</sup> )	エンドキャッピング	pH 範囲
ODS-3	3, 5, 10	100	1.0	450	15.5, Monomeric	—	Yes	2.0-9.0
ODS-2	5	150	1.1	310	18.5, Monomeric	3.50	Yes	2.0-9.0
C8	5	150	1.1	310	12.6, Monomeric	5.00	Yes	2.0-9.0
Phenyl	5	100	—	450	10.0, Polymeric	—	No	2.0-9.0
Silica	5	100	—	450	—	—	No	—

## ■ Prodigy ODS-3

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
3	2.0	50	00B-4222-B0	60,000	AJ0-4286	10	43,000
3	2.0	100	00D-4222-B0	65,000	AJ0-4286	10	43,000
3	2.0	150	00F-4222-B0	70,000	AJ0-4286	10	43,000
3	3.0	100	00D-4222-Y0	65,000	AJ0-4286	10	43,000
3	3.0	150	00F-4222-Y0	70,000	AJ0-4286	10	43,000
3	4.0	100	00D-4222-D0	65,000	AJ0-4287	10	43,000
3	4.6	30	00A-4222-E0	40,000	AJ0-4287	10	43,000
3	4.6	50	00B-4222-E0	50,000	AJ0-4287	10	43,000
3	4.6	100	00D-4222-E0	60,000	AJ0-4287	10	43,000
3	4.6	150	00F-4222-E0	65,000	AJ0-4287	10	43,000
3	4.6	250	00G-4222-E0	67,000	AJ0-4287	10	43,000
5	2.0	50	00B-4097-B0	50,000	AJ0-4286	10	43,000
5	2.0	150	00F-4097-B0	70,000	AJ0-4286	10	43,000
5	2.0	250	00G-4097-B0	75,000	AJ0-4286	10	43,000
5	3.0	150	00F-4097-Y0	70,000	AJ0-4286	10	43,000
5	3.0	250	00G-4097-Y0	72,000	AJ0-4286	10	43,000
5	3.2	150	00F-4097-R0	70,000	AJ0-4287	10	43,000
5	3.2	250	00G-4097-R0	75,000	AJ0-4287	10	43,000
5	4.6	30	00A-4097-E0	40,000	AJ0-4287	10	43,000
5	4.6	50	00B-4097-E0	50,000	AJ0-4287	10	43,000
5	4.6	100	00D-4097-E0	55,000	AJ0-4287	10	43,000
5	4.6	125	00E-4097-E0	58,000	AJ0-4287	10	43,000
5	4.6	150	00F-4097-E0	60,000	AJ0-4287	10	43,000
5	4.6	250	00G-4097-E0	62,000	AJ0-4287	10	43,000
5	10.0	250	00G-4097-N0	200,000	AJ0-7221	3	58,000
10	3.9	300	00H-4244-C0	80,000	AJ0-4287	10	43,000
10	4.6	250	00G-4244-E0	62,000	AJ0-4287	10	43,000
10	10.0	250	00G-4244-N0	200,000	AJ0-7221	3	58,000

## ■ Prodigy ODS-3V (バリデーションサポート)

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
3	4.6	100	00D-4243-E0	60,000	AJ0-4287	10	43,000
3	4.6	150	00F-4243-E0	65,000	AJ0-4287	10	43,000
5	4.6	150	00F-4241-E0	65,000	AJ0-4287	10	43,000
5	4.6	250	00G-4241-E0	67,000	AJ0-4287	10	43,000
10	4.6	250	00G-4088-E0	お問い合わせ	AJ0-4287	10	43,000

## ■ Prodigy ODS-2

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
5	2.0	150	00F-3300-B0	70,000	AJ0-4286	10	43,000
5	3.0	250	00G-3300-Y0	72,000	AJ0-4286	10	43,000
5	3.2	150	00F-3300-R0	70,000	AJ0-4287	10	43,000
5	3.2	250	00G-3300-R0	75,000	AJ0-4287	10	43,000
5	4.6	30	00A-3300-E0	40,000	AJ0-4287	10	43,000
5	4.6	50	00B-3300-E0	50,000	AJ0-4287	10	43,000
5	4.6	75	00C-3300-E0	53,000	AJ0-4287	10	43,000
5	4.6	100	00D-3300-E0	55,000	AJ0-4287	10	43,000
5	4.6	150	00F-3300-E0	60,000	AJ0-4287	10	43,000
5	4.6	250	00G-3300-E0	62,000	AJ0-4287	10	43,000

## ■ Prodigy C8

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
5	2.0	50	00B-3301-B0	50,000	AJ0-4289	10	43,000
5	2.0	150	00F-3301-B0	70,000	AJ0-4289	10	43,000
5	3.0	150	00F-3301-Y0	70,000	AJ0-4289	10	43,000
5	3.0	250	00G-3301-Y0	72,000	AJ0-4289	10	43,000
5	4.0	125	00E-3301-D0	55,000	AJ0-4290	10	43,000
5	4.6	30	00A-3301-E0	40,000	AJ0-4290	10	43,000
5	4.6	50	00B-3301-E0	50,000	AJ0-4290	10	43,000
5	4.6	100	00D-3301-E0	55,000	AJ0-4290	10	43,000
5	4.6	150	00F-3301-E0	60,000	AJ0-4290	10	43,000
5	4.6	250	00G-3301-E0	62,000	AJ0-4290	10	43,000

## ■ Prodigy Phenyl (PH-3)

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
5	4.6	100	00D-4298-E0	55,000	AJ0-4351	10	43,000
5	4.6	150	00F-4298-E0	60,000	AJ0-4351	10	43,000
5	4.6	250	00G-4298-E0	62,000	AJ0-4351	10	43,000

## ■ Prodigy Silica

粒子径 (μm)	内径 (mm)	長さ (mm)	P/N	希望小売価格 (税抜 / 円)	ガードカートリッジ		
					P/N	入数	希望小売価格 (税抜 / 円)
5	4.6	100	00D-4098-E0	55,000	AJ0-4348	10	43,000
5	4.6	250	00G-4098-E0	62,000	AJ0-4348	10	43,000
10	4.6	250	00G-4245-E0	62,000	AJ0-4348	10	43,000

ガードカートリッジをご使用の際は P.87 掲載の専用ホルダーをお求めください。