

(別紙)
 代表的な製品を掲載しております。
 未掲載製品につきましてはお問い合わせください。

2023年12月11日
 株式会社島津ジーエルシー
 第M23-42-02号

取扱中止となる VICI Jour社製品一覧
実施日: 即日

販売中止品		代替品
P/N (製品番号)	製品名	P/N (製品番号)
JR-5FR1TI-5	Frit, Ti, 1/16" diam. x 0.062" thickn., 0.5 um, 5/pkg	(お問い合わせください)
JR-10FR8TI-5	Frit, Ti, 1/2" diam. x 0.040" thickn., 10 um, 5/pkg	(お問い合わせください)
JR-2FR1KTI-5	Frit, Ti, 1" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-2FR1TI-5	Frit, Ti, 1/16" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-2FR2TI-5	Frit, Ti, 1/8" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-2FR4TI-5	Frit, Ti, 1/4" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-2FR6TI-5	Frit, Ti, 3/8" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-2FR8TI-5	Frit, Ti, 1/2" diam. x 0.040" thickn., 2 um, 5/pkg	(お問い合わせください)
JR-C360C	Cap, 360um, SS	C360C
JR-C360CFS	Cap, 360um, SS, for FS tubing > 10'000 psi	C360CFS
JR-C360ET	Make up extender tool, for 360um nuts	C360ET
JR-C360NFFS	Nut, 360um, No Twist, for FS tubing (for > 10 kpsi)	C360NFFS
JR-C360NFPK	Nut, 360um, No Twist, PEEK, with PEEK ferrule	(お問い合わせください)
JR-C360NFS6	Nut, 360um, No Twist, SS, with 316 ferrule	C360NFS6
JR-C360PPK	Plug, 360um, PEEK	C360PPK
JR-C360QTPK4	Tee, 360um, quick mount, 100um bore, for PEEK and FS Tubing (for < 10 kpsi)	(お問い合わせください)
JR-C360QTPK6	Tee, 360um, quick mount, 150um bore, for PEEK and FS Tubing (for < 10 kpsi)	(お問い合わせください)
JR-C360RU.5FS2	Reducing union 1/32" to 360um, 50um bore, for FS tubing (for > 10 kpsi)	C360RU.5FS2
JR-C360RU.5FS4	Reducing union 1/32" to 360um, 100um bore, for FS tubing (for > 10 kpsi)	C360RU.5FS4
JR-C360RU.5FS6	Reducing union 1/32" to 360um, 150um bore, for FS tubing (for > 10 kpsi)	C360RU.5FS6
JR-C360RU.5S62	Reducing union 1/32" to 360um, 50um bore, for metal tubing (up to 20 kpsi)	C360RU.5S62
JR-C360RU.5S64	Reducing union 1/32" to 360um, 100um bore, for metal tubing (up to 20 kpsi)	C360RU.5S64
JR-C360RU.5S66	Reducing union 1/32" to 360um, 150um bore, for metal tubing (up to 20 kpsi)	C360RU.5S66
JR-C360RU1PK6	Reducing union 1/16" to 360um, 150um bore, for PEEK+FS tubing (for < 10 kpsi)	C360RU1PK6
JR-C360RU1S66	Reducing union 1/16" to 360um, 150um bore, for metal tubing (up to 20 kpsi)	C360RU1S66

(別紙)
 代表的な製品を掲載しております。
 未掲載製品につきましてはお問い合わせください。

2023年12月11日
 株式会社島津ジーエルシー
 第M23-42-02号

JR-C360UFS4	Union, 360um, 100um bore, for FS tubing (for > 10 kpsi)	C360UFS4
JR-C360UPK2	Union, 360um, 50um bore, for PEEK and FS tubing (for < 10 kpsi)	(お問い合わせください)
JR-C360UPK4	Union, 360um, 100um bore, for PEEK and FS tubing (for < 10 kpsi)	(お問い合わせください)
JR-C360UPK6	Union, 360um, 150um bore, for PEEK and FS tubing (for < 10 kpsi)	(お問い合わせください)
JR-C-EN.5FPKB	Nut, PEEK, External, 1/32"	C-EN.5FPKB
JR-CFE-S10-5	In-line filter Cartridge, 10 um, 5/pkg	CFE-S10
JR-CFE-S2-5	In-line filter Cartridge, 2 um, 5/pkg	CFE-S2
JR-CFE-S75-5	In-line filter Cartridge, 75 um, 5/pkg	CFE-S75
JR-C-NERU1FPK	Union, PEEK, 150 um bore, reducing, Ext.1/16" - Int. 1/32"	C-NERU1FPK
JR-C-NEU.5FPK	Union, PEEK, 150 um bore, Int. 1/32" - Ext. 1/32" ZDVF	C-NEU.5FPK
JR-C-NEU.5XFPK	Union, PEEK, 100 um bore, Int. 1/32" - Ext. 1/32" ZDVF	C-NEU.5XFPK
JR-C-NNFFPK	Fitting, PEEK, one-piece short, grooved ferrule, 1/32"	C-NNFFPK
JR-EN1FPKB	Nut, PEEK, External, 1/16"	(お問い合わせください)
JR-FRK1	Ferrule removal kit	FRK1
JR-LZN1-5	Nut, SS, 1/16" Long, Valco type, 10-32, 5/pkg	LZN1-10
JR-MZN1-5	Nut, SS, 1/16" Medium, Valco type, 10-32, 5/pkg	MZN1-10
JR-ZC.5FPK	Cap 1/32", fingertight, PEEK	ZC.5FPK
JR-ZC1	Cap, SS, female 1/16"	ZC1
JR-ZC1PK	Cap, PEEK, female 10-32	(お問い合わせください)
JR-ZC2	Cap, SS, female 1/8"	ZC2
JR-ZF.5PK-5	Ferrule, PEEK, 1/32", 5 /pkg	ZF.5PK-5
JR-ZF.5S6-5	Ferrule, 316SS, 1/32", Valco type, 5/pkg	ZF.5S6-10
JR-ZF1S6-5	Ferrule, 316SS, 1/16", Valco type, 5/pkg	ZF1S6-10
JR-ZF2PK-5	Ferrule, PEEK, 1/8", 5/pkg	ZF2PK-10
JR-ZF2S6-5	Ferrule, 316SS, 1/8", Valco type, 5/pkg	ZF2S6-10
JR-ZN.5-5	Nut, SS, 1/32", Valco type, 6-40, 5/pkg	ZN.5-10
JR-ZN.5FPK-5	Nut, PEEK, fingertight, 1/32", 5 /pkg	ZN.5FPK-10
JR-ZN1-5	Nut, SS, 1/16", Valco type, 10-32, 5/pkg	ZN1-10

(別紙)
 代表的な製品を掲載しております。
 未掲載製品につきましてはお問い合わせください。

2023年12月11日
 株式会社島津ジーエルシー
 第M23-42-02号

JR-ZN2-5	Nut, SS, 1/8", Valco type, 5/16"-24, 5/pkg	ZN2-10
JR-ZN2PK	Nut, PEEK, Hex-Head 1/8"	ZN2PK-10
JR-ZN2PK-5	Nut, PEEK, Hex-Head 1/8", 5/16"-24, 5/pkg	ZN2PK-10
JR-ZP.5FPK	Plug 1/32", fingertight, PEEK	ZP.5FPK
JR-ZP1FPK	Plug 1/16", fingertight, PEEK	ZP1FPK
JR-ZRU1.5FPK	Union, PEEK, reducing 1/16" - 1/32", 10-32, bore 0.25mm, complete	ZRU1.5FPK
JR-ZRU21	Union, SS, reducing 1/8" to 1/16", bore 0.75mm, complete	ZRU21
JR-ZRU21FPK	Union, PEEK, HP reducing complete, 1/8" to 1/16", thru-hole 0.75 mm	ZRU21FPK
JR-ZT1	Tee, SS, 0.75 mm bore, 1/16", 10-32, complete	ZT1
JR-ZT1C	Tee, SS, 0.25 mm bore, 1/16", 10-32, complete	ZT1C
JR-ZT1M	Tee, SS, 0.5 mm bore, 1/16", 10-32, complete	ZT1M
JR-ZU1CS6	Union, SS, bore 0.25 mm, 1/16", 10-32, complete	(お問い合わせください)
JR-ZU1MS6	Union, SS, bore 0.50mm, 1/16", 10-32, complete	(お問い合わせください)
JR-ZU1S6	Union, SS, bore 0.75mm, 1/16", 10-32, complete	(お問い合わせください)
JR-ZU1XCS6	Union, SS, bore 0.15 mm, 1/16", 10-32, complete	(お問い合わせください)
JR-ZX1	Cross, SS, 0.75 mm bore, 1/16", 10-32, complete	ZX1
JR-ZX1C	Cross, SS, 0.25mm bore, 1/16", 10-32, complete	ZX1C
JR-ZX1M	Cross, SS, 0.50mm bore, 1/16", 10-32, complete	ZX1M